

ELEA SIT-TO-STAND AID – TRANSFER AID

Made by Leonair



- Enables healthy-back working
- Easy to operate
- Provides motorised support
- Highly tilt-resistant
- Avoids shoulder problems due to belt attachment in pelvic area
- Promotes existing mobility, also suitable for standing exercises
- Provides firm and safe support during transfer



Basic equipment

- Support column with a platform
- 2 wheels with a locking brake
- 4 wheels without a locking brake
- Multi-grip bars
- Shinbone rest, adjustable in height and width
- Patient belt
- Lift system (motuna) providing motorised support, including manual control unit, battery pack and charger



CE






HMV: 22.40.01.1022

Leonair®


Indication	Persons with a limited range of motion and functional restrictions in the lower extremities and the pelvic girdle, which make it difficult or impossible for them to get up without assistance.	
Contraindication	Persons unable to put weight on their feet for a short time although having handlebars for both hands and being supported with a tension band and shinbone rests.	
Specifications	Electric pulling system (System motuna)	24 VDC / 2.9 A Safety nut in push direction, brake in push direction
	Battery module	Battery pack consisting of 2 rechargeable batteries
	Wall charging unit	externally connectable / suitable for the aforementioned battery pack
	Total height	130 cm (without accessories)
	Total length	85 cm
	Construction width	60 cm
	Ground clearance	8 cm
	Front wheels	2 x 100 mm with a locking brake
	Back wheels	4 x 50 mm
	Weight	45 kg
	Loading capacity	130 kg

Order data	Order no.
Basic version ELEA	09.100.001

Zubehör

	Patient belt, size 1, maximal abdominal circumference 100 cm	00.011.926
	Patient belt, size 2, maximal abdominal circumference 130 cm	00.011.927
	Chest rest. Allows patients to gently rest their chest or chin, e.g. while getting cleaned after using the toilet.	00.012.260
	Vertical rest including chest rest (e.g. for standing exercises). Supports upright position.	00.012.261
	Foot holders including instep straps, maximal shoe size 37 (pair). Attaches feet to standing surface. May be used for patients with spasms.	00.012.262
	Foot holders including instep straps, maximal shoe size 45 (pair)	00.012.263
	Calf holders (pair) To stabilise legs and knees if these tend to give way laterally.	00.011.860
	Additional handles (pair) To provide more options for patients to hold on to.	00.012.264
	Padded inserts for shinbone rests (pair)	00.012.436

The accessories can also be used for the ROLEO stand-up aid.



Further accessories and customisation are possible upon request. Our ELEA sit-to-stand and transfer aid meets the requirements according to the European Medical Device Regulation 2017/745 „EU MDR“. We reserve the right to make technical changes and changes in colour at any time.

HMV: 22.40.01.1022