



Enables healthy-back working – easy to operate – avoids shoulder problems due to belt attachment in the body's centre of gravity

- Electric transfer aid
- Supports the upper body
- Handles may be telescoped and pivoted
- Adjustable patient belt
- Adjustable shinbone rests
- Wheeled aid

TRANSFER AID

For care-dependent people with or without core stability

| | | |
|-------------------|---|---|
| Indication | Persons with a limited range of motion and functional restrictions, which make it difficult or impossible for them to get up without assistance. | |
| Contraindication | Persons unable to put weight on their legs and/or feet for a short time although their upper body is supported and their legs are fixed with shin-bone rests. | |
| Specifications | Electric pulling system | 24 VDC / 2,9 A Safety nut in push direction, brake in push direction |
| | Battery module | Battery pack consisting of 2 rechargeable batteries |
| | Wall charging unit | externally connectable |
| | Total height | 112 cm |
| | Total length | 92 cm / handles retracted |
| | Construction width | 59 cm |
| | Ground clearance | 6.5 cm |
| | Back wheels | 2 x Ø 100 mm with a locking brake |
| | Front wheels | 3 x Ø 35 mm |
| | Weight | 33 kg |
| | Loading capacity | 130 kg |
| Order data | | Order no. |

| | |
|--|------------|
| Basic version | 10.100.001 |
| Patient belt, size 1, maximal abdominal circumference 100 cm | 00.012.428 |
| Patient belt, size 2, maximal abdominal circumference 130 cm | 00.012.429 |
| Padded inserts for shinbone rests (pair) | 00.012.436 |
| Pulling belt | 00.012.469 |

Further accessories upon request: foot holders, calf holders; customisation is possible.

Our LiLLY transfer aid meets the requirements according to the European Medical Device Regulation 2017/745 „EU MDR“. We reserve the right to make technical changes and changes in colour at any time.




Abb. zeigt Zubehör