

ROLEO® SIT-TO-STAND AID – TRANSFER AID

Made by Leonair



Turning around with ROLEO makes transfer easy and saves space.

- Enables healthy-back working
- Easy to operate
- Provides motorised support
- Highly tilt-resistant
- Easy on joints
- Avoids shoulder problems due to belt attachment in pelvic area
- Promotes existing mobility, also suitable for standing exercises
- Provides firm and safe support during transfer









Basic equipment

- Support column
- Anti-tilt supports at front and back
- 4 wheels, locking brake
- 360° turntable
- Multi-grip bars
- Shinbone rest, adjustable
- Patient belt
- Lift system providing motorised support, including manual control unit, battery pack and charger



HMV: 22.40.01.1023

Leonair®

Indication	Persons with a limited range of motion and functional restrictions in the lower extremities and the pelvic girdle, which make it difficult or impossible for them to get up without assistance. Patients confined to bed, for short-term seated transfer and transport in sitting position.	
Contraindication	Persons unable to put weight on their feet for a short time although having handlebars for both hands and being supported with a tension band and shinbone rests.	
Specifications	Electric pulling system (System motuna)	24 VDC / 2,9 A Safety nut in push direction, brake in push direction
	Battery module	Battery pack
	Charger	suitable for the aforementioned battery pack
	Working area turntable	360°, screened in 32 steps
	Total height	127 cm (without accessories)
	Total length	93 cm standard / with a seat 105 cm
	Construction width	59 cm
	Ground clearance	14 cm
	Wheels	4 x 100 mm with a locking brake
	Weight	54 kg
	Loading capacity	130 kg
Order data		Order no.
     	Basic version ROLEO	08.502.002
	Patient belt, size 1, max. abdominal circumference 100 cm	00.011.926
	Patient belt, size 2, max. abdominal circumference 130 cm	00.011.927
	Chest rest	00.012.260
	Vertical rest including chest rest	00.012.261
	Foot holders incl. instep straps, max. shoe size 37 (pair)	00.012.262
	Foot holders incl. instep straps, max. shoe size 45 (pair)	00.012.263
	Calf holders (pair)	00.011.860
	Additional handles (pair)	00.012.264
	Seat for safe and comfortable patient transport, e.g. bed to toilet without a wheelchair	00.012.265
	Padded inserts for shinbone rests (pair)	00.012.436

The accessories, except for the seat, can also be used for the ELEA stand-up aid.

Further accessories and customisation are possible upon request. Our ROLEO sit-to-stand and transfer aid meets the requirements according to the European Medical Device Regulation 2017/745 „EU MDR“. We reserve the right to make technical changes and changes in colour at any time.

HMV: 22.40.01.1023